


# AM-U-SeatSling : Universal Seated Transfer


## TRIM INDICATES SIZE/WEIGHT\*

<b>S</b>	31-75 LBS	14-34 KG
<b>M</b>	75-145 LBS	34-66 KG
<b>L</b>	145-220 LBS	66-100 KG
<b>XL</b>	220-440 LBS	100-200 KG
<b>XXL</b>	440-800 LBS	200-363 KG

\* IF IN DOUBT STEP DOWN ONE SIZE




**1** **TO PLACE SLING UNDER PATIENT:** ELEVATE HEAD OF THE BED. LEAN PATIENT FORWARD. PLACE SLING BEHIND PATIENT WITH LABEL ON THE OUTSIDE. THE LEG WINGS SHOULD COME ALONG THE SIDE OF THE PATIENT.



**2** USE BACK POCKET TO SLIDE CENTER OF SLING DOWN TO THE COCCYX AREA



**3** FOLD EACH LEG WING IN TWO AND PUSH UNDERNEATH PATIENT'S LEGS JUST ABOVE THE KNEE.



**4** TAKE LEFT WING AND PULL IT THROUGH DIGNITY STRAP LOCATED ON THE RIGHT LEG.


## AM-U-SeatSling : Universal Seated Transfer (cont.)



5


**TO LIFT USING THE CLIPS:** POSITION HANGER BAR IN FRONT OF PATIENT. PLACE LOWER PART OF CLIP HOLE OVER HANGER BAR KNOB, SLIDE CLIP DOWN UNTIL IT LOCKS. REPEAT ON OPPOSITE SIDE

6




AFTER LIFTING, TRANSPORT PATIENT BY HOLDING FRONT HANDLE ON HANGER BAR. PLACE PATIENT

7



BEGIN TO LOWER PATIENT AND USE FRONT HANDLE TO PUSH INTO CHAIR BAR.



5b

**TO LIFT USING THE LOOPS:** ATTACH STRAPS TO HANGER BAR  
TUCK CLIPS IN THEIR POCKETS.

FOR SEATED POSITION: SHORT LOOP FOR CHEST STRAPS AND LONG LOOP FOR LEG STRAPS

FOR SEMI-RECLINED TRANSFER: LONG LOOP FOR CHEST STRAPS AND SHORT LOOP FOR LEG STRAPS

6b

AFTER LIFTING, MOVE PATIENT BY HOLDING THE SLING OR THE HANDLES LOCATED BEHIND THE SLING.

7b

TO LOWER PATIENT INTO CHAIR, PUSH GENTLY FORWARD ON THE LEG WINGS TO ENSURE PATIENT'S BUTTOCKS ARE IN CONTACT WITH THE BACK OF THE SEAT.